

Back Support Order / Quote Request

March 2013

circle one: **NEW ORDER** **QUOTE REQUEST**

Company Name: _____	PO #: _____	Mark For: _____
Contact Person: _____	Phone #: _____	Send Order Confirmation To: _____
Delivery Address: _____		

Fax to 210.341.0009 or email to cs@adirides.com

Build your back by circling applicable specs...

<u>Aluminum Shell</u>	<u>Carbon Shell</u>	<u>Back Height</u>	<u>Contour</u>	<u>Hardware</u>
10" / 25.4cm	13" / 33cm	Low 10" / 25.4cm	Active 2.5" / 63.5mm	Fixed Elite (carbon only)
11" / 28cm	14" / 35.5cm	Med 13" / 33cm	Deep 4" / 100mm	2 Point PRO
12" / 30.5cm	15" / 38cm	Std 16" / 40.6cm	Extra Deep 6" / 100mm (aluminum only)	Quick Release
13" / 33cm	16" / 40.6cm	Tall 18" / 45.7cm (aluminum only)	Other _____	Width Adj. Elite
14" / 35.5cm	17" / 43cm	Extra Tall 20" / 51cm (aluminum only)	_____	
15" / 38cm	18" / 45.7cm	Other _____		
16" / 40.6cm	19" / 48cm			
17" / 43cm				
18" / 45.7cm				
19" / 48cm				
20" / 51cm				
Other _____				

<u>Lateral Support Options</u>	
<u>Pad Size</u>	<u>Hardware</u>
Small 3" H x 3" D	Stationary
Med 3" H x 4" D	Swing-Out
Large 4.5" H x 4.5" D	
X Large 4" H x 6" D	
Other _____	

<u>Back Support Accessories</u>
Lumbar Support
Privacy Flap
Loop Handle (quick release only)

PRODUCT ITEM NUMBER FOR THIS REQUEST - _____

ITEM NUMBER BUILDER - All item number begin with shell type: AL (for Aluminum) or CF (Carbon Fiber) + BR

To build the correct product item number for your request use the following formula:

AL or CF + BR + Wheelchair Width in inches + Back Support Height (L, M, S, T or X) + Contour (A, D or 6) + Hardware (F, 2, R or W)

Eg: An Aluminum Series Back Support to fit a 10" wide wheelchair, back height low with active contour and quick release mounting has the item number **ALBR10LAR**



401 Isom Road, Suite 520, San Antonio, Texas 78216
 p: 210.341.0008 f: 210.341.0009 t: 888.684.2234
 e: cs@adirides.com w: adirides.com