




# Back Support Shell Specifications

## Contour Width & Depth

Active Contour																						
width	10"		11"		12"		13"		14"		15"		16"		17"		18"		19"		20"	
	inch	cm	inch	cm	inch	cm	inch	cm	inch	cm	inch	cm	inch	cm	inch	cm	inch	cm	inch	cm	inch	cm
<b>A</b>	7.5	19.1	8.5	21.6	9.5	24.1	10.5	26.7	11.375	28.9	12.25	31.1	13.5	34.3	14.375	36.5	15.375	39.1	16.375	41.6	17.25	43.8
<b>B</b>	7	17.8	8	20.3	9	22.9	10	25.4	10.75	27.3	11.75	29.8	12.75	32.4	13.75	34.9	14.75	37.5	15.75	40.0	16.75	42.5
<b>C</b>	2.1	5.3	2.1	5.3	2.1	5.3	2.1	5.3	2.25	5.7	2.25	5.7	2.25	5.7	2.25	5.7	2.25	5.7	2.25	5.7	2.25	5.7

Deep Contour																						
width	10"		11"		12"		13"		14"		15"		16"		17"		18"		19"		20"	
	inch	cm	inch	cm	inch	cm	inch	cm	inch	cm	inch	cm	inch	cm	inch	cm	inch	cm	inch	cm	inch	cm
<b>A</b>	7.5	19.1	8.5	21.6	9.5	24.1	10.5	26.7	12.625	32.1	13.75	34.9	14.25	36.2	15.25	38.7	16.75	42.5	17.75	45.1	18.5	47.0
<b>B</b>	6	15.2	7	17.8	8	20.3	9	22.9	11	27.9	12	30.5	13	33.0	14	35.6	15	38.1	16	40.6	17	43.2
<b>C</b>	3.625	9.2	3.625	9.2	3.625	9.2	3.625	9.2	4.25	10.8	4.25	10.8	4.25	10.8	4.25	10.8	4.25	10.8	4.25	10.8	4.25	10.8

Extra Deep Contour										
width	16"		17"		18"		19"		20"	
	inch	cm	inch	cm	inch	cm	inch	cm	inch	cm
<b>A</b>	15.5	39.4	16.5	41.9	17.5	44.5	18.5	47.0	19.5	49.5
<b>B</b>	14.75	37.5	15.75	40.0	16.75	42.5	16.75	42.5	17.75	45.1
<b>C</b>	6	15.2	6	15.2	6	15.2	6	15.2	6	15.2



- A. Actual Back Width
- B. Trunk Width
- C. Depth