

## ZIP Privacy Flap

ADI privacy flaps are an optional feature. Our new design simply “Zips” onto the bottom of the fabric cover for a secure hold. This design allows you to add a privacy flap at the time of order or when needed.

### STEP #1

attach privacy flap to the fabric back support cover with zipper



### STEP #2

Velcro privacy flap to seatbase



Note:

a strip of Velcro hook & loop is sewn into flap to attach to all bases / slings



Fabric Cover

Privacy Flap Zipper

ZIP Privacy Flap

