

Low Deep (10" height / 4" contour)											
Width	10"	11"	12"	13"	14"	15"	16"	17"	18"	19"	20"
A	4.25	4.5	5	5.25	5.5	5.5	5.5	5.5	5.5	5.5	5.5
B	5	5	6.25	6.25	6.25	6.25	6.25	6.25	6.25	6.25	6.25
C	10	10	10	10	10	10	10	10	10	10	10

Medium Deep (13" height / 4" contour)									
Width	12"	13"	14"	15"	16"	17"	18"	19"	20"
A	6.5	6.5	6.75	6.75	6.75	6.75	6.75	6.75	6.75
B	7	7	7	7	7	7	7	7	7
C	13	13	13	13	13	13	13	13	13

Standard Deep (16" height / 4" contour)							
Width	14"	15"	16"	17"	18"	19"	20"
A	7.25	7.25	7.25	7.25	7.25	7.25	7.25
B	7	7	7	7	7	7	7
C	16	16	16	16	16	16	16

Standard Extra Deep (16" height / 4" contour)					
Width	16"	17"	18"	19"	20"
A	7.25	7.25	7.25	7.25	7.25
B	6	6	6	6	6
C	16	16	16	16	16

Tall Deep (18" height / 4" contour)					
Width	16"	17"	18"	19"	20"
A	9.5	9.5	9.5	9.5	9.5
B	7	7	7	7	7
C	18	18	18	18	18

Tall Extra Deep (18" height / 6" contour)					
Width	16"	17"	18"	19"	20"
A	9.5	9.5	9.5	9.5	9.5
B	6	6	6	6	6
C	18	18	18	18	18

Extra Tall Deep (20" height / 6" contour)					
Width	16"	17"	18"	19"	20"
A	11	11	11	11	11
B	7	7	7	7	7
C	20	20	20	20	20

Extra Tall Extra Deep (20" height / 6" contour)					
Width	16"	17"	18"	19"	20"
A	11	11	11	11	11
B	6	6	6	6	6
C	20	20	20	20	20



Back Support Deep Contour Height Specs

Note:
All dimensions are in inches

