



# Backrest Basics



item #:

**A** **L** **B** **R** **1** **0** **L** **A** **R**

shell type      back rest      chair width      height      contour      mounting

AL - AL Series

CF - CF Series

GU - G Series

note:

Shell Assembly Weight

AL Series - 1.5lb / 680gm

CF Series - .820lb / 370gm

- above weights are of a 16" STD

- all weights approx.

10" / 25cm

11" / 28cm

12" / 30.5cm

13" / 33cm

14" / 35.5cm

15" / 38cm

16" / 40.5cm

17" / 43cm

18" / 45.5cm

19" / 48cm

20" / 51cm

L

low  
(10"/25cm)

M

medium  
(13"/33cm)

S

standard  
(16"/40.5cm)

T

tall  
(18"/45.7cm)

X

x-tall  
(20"/51cm)

A

active  
(2.5"/63.5mm)

D

deep  
(4"/100mm)

6

extra  
deep  
(6"/150mm)

F

fixed elite  
(.60lb / 272gm)

2

2-point pro  
(1.25lb/566gm)

R

quick  
release  
(1.30lb / 589gm)

W

width  
adjustable  
elite  
(.63lb / 285gm)

E

eq release  
(1.50lb / 680gm)

Accessible Designs, Inc.  
401 Isom Road, Suite 520  
San Antonio, Texas 78216

e: cs@adirides.com ~ w: adirides.com

