




Back Support Contour Specs

Measurements on Centimeters

Active Contour											
Chair width	25.40	27.94	30.48	33.02	35.56	38.10	40.64	43.18	45.72	48.26	50.80
A	20.96	23.50	26.04	28.57	33.02	35.56	38.10	40.64	43.18	45.72	48.86
B	17.78	20.32	22.86	25.40	28.57	31.12	33.66	36.20	38.73	41.27	43.81
C	6.35	6.35	6.35	6.35	6.35	6.35	6.35	6.35	6.35	6.35	6.35

Deep Contour											
Chair width	25.40	27.94	30.48	33.02	35.56	38.10	40.64	43.18	45.72	48.26	50.80
A	20.96	23.50	26.04	28.57	34.92	37.47	40.01	42.55	45.09	47.63	50.16
B	17.14	19.68	22.23	24.77	29.21	31.75	34.29	36.83	39.37	41.91	44.45
C	10.16	10.16	10.16	10.16	10.16	10.16	10.16	10.16	10.16	10.16	10.16

Extra Deep Contour					
Chair width	40.64	43.18	45.72	48.26	50.80
A	40.01	42.55	45.09	47.63	50.16
B	32.38	34.92	37.47	40.01	42.55
C	15.24	15.24	15.24	15.24	15.24



- A. Actual Back Width
- B. Trunk Width
- C. Depth