


For use with: Fusion16 DB Series Wheels

NOTE: for the above wheels, axle needs to be adjusted to 3.875" / 98.5mm in length as shown below




Item #: AXL005

Description: ADI quick release axle, length 96.5-103mm

(USE: Fusion DB, Fusion Everyday and Spinerger)

For use with: Spinerger Spox & LX (24" & 25") w/ Spline Drive Insert

NOTE: for the above wheels, axle needs to be adjusted to 4.100" / 104mm in length as shown below



Item #: AXL004

Description: ADI quick release axle, length 103-109.5mm

(USE: Spinerger w/ DB spline drive insert only)